<i>Osobní váha s tělesnou analýzou • NÁVOD K OBSLUZE</i>	<b>CZ</b>	4-12
<i>Osobná váha s telesnou analýzou • NÁVOD NA OBSLUHU</i>	<b>SK</b>	13-21
<i>Body analysis scale • USER MANUAL</i>	<b>EN</b>	22-29
<i>Személyi mérleg test analízissal • HASZNÁLATI ÚTMUTATÓ</i>	<b>HU</b>	30-37
<i>Waga łazienkowa z analizą ciała • INSTRUKCJA OBSŁUGI</i>	<b>PL</b>	38-45
<i>Personenwaage mit Körperanalyse • BEDIENUNGSANLEITUNG</i>	<b>DE</b>	46-54

# LENA


1

2

3

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HOUSEHOLD USE ONLY. DO NOT IMMERSE IN WATER OR OTHER LIQUIDS.

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UNIT - jednotka; OPEN - otvorit

**TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, ARRANGES OR PLAYPENS. THIS BAG IS NOT A TOY.**

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
LENA

## INSTRUCTIONS FOR USE

Dear customer, thank you for purchasing our product. Keep this instruction for use and cash receipt (if possible, including the package and its interior) for further reference.

## I. SAFETY WARNING

- Consider the instructions for use as a part of the appliance and pass them on to any other user of the appliance.
- **The appliance is intended for use in households only and for similar use! It is not intended for use in health-care facilities and for commercial use!**
- Never use the appliance if it malfunctions, if it fell to the ground and got damaged or if it fell into water. In these cases bring the appliance to a specialised service to have the safety and functions of the device checked.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- **We recommend: Do not step onto the scale with wet feet or on wet surface of the scale, there is a risk of slipping.**
- **Never immerse the scale into water (even their parts)!**
- If the scale was stored in lower temperatures, let it adjust to the room temperature.
- Use the scale in a horizontal position only and at places where they cannot knock over, in sufficient distance from heat sources (e.g. **stove, fireplace, heat radiators**), sources of increased air humidity (e.g. **sauna, bathroom, swimming pool**) and equipment with strong electromagnetic field (e.g. **microwave oven, mobile phone**).
- Protect the personal scale from dust, chemicals, excessive humidity, toiletries, liquid cosmetic preparations, intensive temperature changes and do not expose them to direct sunlight.
- Do not step on the edge of the scale! The scale could overturn and cause injury!
- Do not place the scale on uneven surface or on surface covered with carpet. Unstable or soft surface under the scale can have negative effect on precision of weighing.
- In handling the scale, make sure injury of persons or damage of the appliance is prevented.
- Do not overload the appliance with weight higher than the maximum capacity! You will thus prevent damaging the appliance
- Do not use the fat analysis feature on this scale if you are using one of the following devices: pace maker, electrocardiograph and/or other medical device implanted in the body or device used as support.
- Handle the scale with care. Do not throw them or jump on them. The scale is durable, but improper handling can destroy their electronic sensors.
- Do not put any objects on the surface of the scale during storing, they could be damaged.
- If the battery is leaking, remove it immediately; it can damage the scale.

- Do not use rechargeable accumulators.
- Dispose of a discharged batteries in a suitable way.
- Never use the appliance for any other purpose than for that described in these instructions for use!
- Keep in mind that the data read from this scale are approximate and inaccurate in comparison to results of actual medicinal analyses. The exact determination of biometrical characteristics can be done only by a doctor via medical methods.
- Keep batteries and scale out of reach of children and incapacitated people. A person, who would swallow the batteries, must seek medical help immediately.
- Do not use during pregnancy, with symptom of edema and during dialysis treatment.
- All text in other languages and pictures shown on the packaging, or the product itself, are translated and explained at the end of this language mutation.
- The producer is not responsible for damage caused by improper use of the appliance (e.g. **breaking the board, etc.**) and its guarantee does not cover the appliance in the case of failure to comply with the safety warnings above.

## II. GENERAL INFORMATION

### Body Fat Ratio Chart

Age	Female (%)				Male (%)			
	Low	Healthy	High	Obese	Low	Healthy	High	Obese
<39	<21	21-33	33-39	39>	<8	8-19	19-25	25>
40-59	<23	23-35	35-40	40>	<11	11-22	22-28	28>
60>	<24	24-36	36-42	42>	<13	13-25	25-30	30>

The optimal percentage of body water is about 50-65 % for men and 45-60 % for a women.

### Muscle Ratio Chart

Age	Female %			Male %		
	Low	Healthy	High	Low	Healthy	High
<30	<34	34-43	>43	<42	42-57	>57
31-50	<31	31-38	>38	<40	40-52	>52
51>	<27	27-34	>34	<37	37-48	>48

### BMI

The technology used in these scales, Bio-electrical Impedance Analysis (BIA) is more technologically advanced than the order Body Mass Index (BMI) method where fat % is calculated using only height and weight. BIA technology uses 5 parameters to get a more accurate resolution for the analysis (weight and biological resistance which are automatically measured by the scales together with the three personal parameters entered manually - sex, height and age). This index sets mutual ratio between the body weight and height. It can be used to determine whether the person suffers from underweight, overweight or obesity or whether his/her health condition is normal. Also with regards to heart and vascular diseases, body mass index is a relatively good risk index, particularly in healthy people. However, the index is for reference only. Experts recommend supplementing the calculation by at least measuring the waist, which makes determining disposition of the body fat more precise.

BMI will not provide very reliable information for small children, elderly people and active sportsmen who have a lot of muscle mass as their results can be biased.

BMI	Underfat	Healthy	Overfat	Obese
	< 18.5	18.5-27	28-39	> 40


### Optimum energy consumption KCAL

It represents optimum amount of energy consumption necessary for maintaining the basic vital functions and balanced metabolism. The data on the display show the optimum daily consumption of energy that you give out at your weight, height, age and gender. You should also take in this amount of energy to keep your metabolism balanced. If you intend to reduce your weight, your daily energy intake should be a bit lower. However, remember to keep carbohydrates, fats and proteins balanced in your diet. 239 cal = 1 kJ

### Important notes

- This scale requires you to weigh yourself barefoot for the reason of correct conductivity during the body fat analysis.
- Do not place the scale on an uneven surface or a surface covered with a carpet.
- **Do not use the fat analysis feature on this scale if you are using one of the following devices: pace maker, electrocardiograph and/or other medical device implanted in the body or device used as support.**
- Excessive drinking, eating, exercising, health conditions, used medicaments, women's menstruation cycle etc. influence the results of the analysis.
- Tables and information concerning optimum arrangement of body characteristics is only orientational.

## III. PREPARATION FOR USE

Remove all the packing material and take the scale out. Remove all possible adhesion foils, stick-on labels or paper from the scale. Open the battery cover on the bottom of the scale, insert the battery with the correct polarity and close the cover in the opposite direction. Place the scale on flat, stable and dry surface. Turn on the scale by tapping your foot on the weighing surface and check the weight unit displayed on the display (kg). If you wish to set the weighing unit (kg / st:lb), you can do so by pressing the **"UNIT"** button at the bottom of the scale.

kg = kilogram

1 kg = 2.2 lb

lb = pound (British weight unit)

1 lb = 0.454 kg

st:lb = stone:pound (American/British weight unit)

1 st = 6.356 kg

FT;In = feet;inches

1 FT = 30.48 cm; 1 In = 2.54 cm

### Replacing the batteries

Open the battery cover on the bottom of the scale (pic. 2). Insert a battery with the right polarity and close the cover. Use AAA batteries (4 pcs) for this appliance.

#### Notes

- Do not combine various types of batteries. Do not even combine used and new batteries together.
- Batteries not intended for this purpose must not be charged.
- Remove the batteries if you will not be using the appliance for a long time.
- The supply terminals must not be short-circuited.

## IV. INSTRUCTIONS FOR USE

Your scale is fine electronic equipment. One of the things that are most difficult to weigh is human body as it is always in motion. If you want to obtain reliable results, always try to step with your feet at the same place on the scale and stand still (Fig. 1). If possible, weigh yourself at the same time of the day (best in the morning) and also on the same place to get comparable results. Remember that excessive clothes can increase your weight as well as eating and drinking before weighing. Short-term weight fluctuation is generally attributed to loss of liquids.

### A) Regular weighing mode

Turn the scale on by stepping on the weighing surface. After a moment, the measured weight flashes and the value remains displayed for several seconds. **Then the scale are switched off automatically.**

### B) Weighing and total body analysis mode


**Setting (change) of personal data for the body fat and water analysis.** Before the first analysis, it is necessary to put in personal data (sex, height and age) which will be saved to the memory. It is possible to save data of 10 persons in the scale's memory.

**The input data can be set and gradually changed in the following range:**


Memory	Gender	Height	Age
0-9	male / female	80 – 220 cm 2'07.5" – 7'02.5"	10-80

1. Turn the scale on by stepping on the weighing surface.
2. Press the **"SET"** button.
3. Select the number of the memory (with button **"^"** or **"v"**) to which you want to save the data for the person to be weighed (**P0-P9**).
4. Confirm by pressing the **"SET"** button.
5. Set your sex with the **"^"** or **"v"** button: **"♂"** = male / **"♀"** = female.
6. Confirm by pressing the **"SET"** button.
7. Set your age with the **"^"** or **"v"** button (Age).
8. Confirm by pressing the **"SET"** button.
9. Set your height with the **"^"** or **"v"** button.
10. Confirm by pressing the **"SET"** button.

Now all the personal data is saved (changed). Wait till 0.0 appears. Stand on the scale and wait for results.

First the weight will be considered. After the display, stand still on the scale and then the measurement will begin to analyze. During measurement, the **"o"** symbol appears alternately on the display. The following data will start to show after a few seconds:

In the upper part of the display, the weight and heart rate value will flash alternately. The results will be displayed for a few seconds and then the scale will automatically turn off.

#### Notice

There must not be a contact between feet, calves or thighs. Otherwise the measurement could be inaccurate.

For a repeated analysis of body characteristics (after saving the input data of the person to the memory), proceed as follows:



1. Turn the scale on by stepping on the weighing surface.
2. Press the **"SET"** button.

3. Use the "∧" or "∨" button to select the memory number (P0-P9) under which your settings are stored.
4. Wait for 0.0 to be displayed.
5. Step on the scale and wait for the weigh-in and body analysis.

#### Note

The scale is switched on automatically if you do not step on it within approx. 10 seconds after „0.0“ is displayed. The scales can be turned off also by pressing the button.

### Error messages and their symbols

Symbols displayed	Interpretation
	Battery is dead
	Analysis error
	The measured weight is higher than maximum allowed value of the scale.

## V. TROUBLESHOOTING OF POSSIBLE PROBLEMS

### If you have problems with weight:

- Check whether the batteries are inserted properly.
- Check whether you selected the proper weight unit.
- Check whether the scale is on a flat floor and they do not touch the wall.
- If unusual signs are displayed, remove the insert the battery again.
- If signs appear, such as if the display is flashing or if there is a faulty sign, try to move the scale away from the influence of a source of interference or ensure switching off the source for the time of using the scale.

## VI. STORING

After using the scale, put it into the horizontal position (not vertically), out of reach of children and incapacitated people (see par. I. **SAFETY WARNINGS**).


## VII. MAINTENANCE

Do not use coarse and aggressive cleaning agents (e.g. **sharp objects, scrapers, diluting agents or other solvents**). Clean the board by wiping with a soft wet cloth. Make sure that water does not get inside the appliance.

**More extensive maintenance or maintenance requiring intervention in the inner parts of the appliance must be carried out by professional service!**

**Failure to comply with the instructions of the manufacturer will lead to expiration of the right to repair!**

The information regarding warranty and product repairs you can find on the website [www.eta.cz/servis-eu](http://www.eta.cz/servis-eu).

## VIII. ENVIRONMENTAL PROTECTION


If the dimensions allow, there are labels of materials used for production of packing, components and accessories as well as their recycling on all parts. The symbols specified on the product or in the accompanying documentation mean that the used electrical or electronic products must not be disposed of together with municipal waste. For proper disposal, hand them over at special collection places where they will be accepted free of charge. Suitable disposal of the product can help to maintain valuable natural resources and to prevent possible negative impacts on the environment and human health, which could be the possible consequences of improper waste disposal. Ask for more details at the local authorities or at a collection site. Fines can be imposed for improper disposal of this type of waste in agreement with the national regulations. Remove a discharged battery and dispose of it in a suitable and safe way through the special collection networks. Never dispose of the battery by burning!

## XI. TECHNICAL DATA

Weighing range max	180 kg / 397 lb / 28.3 st
Discrimination:	0.1 kg
Power supply:	4x Batteries 1.5 V AAA (included)
Net weight of device:	2 kg
Memory capacity:	10
Mass alignment:	from 80 to 220 cm (2'07.5" – 7'02.5" „Ft:in)
Age:	from 10 to 100 years
Minimum weight:	from 5 kg (11 lb / 0:8 st)
Maximum weight:	180 kg (397 lb/28.3 st)
Body fat measuring range:	5 % – 70 %
Body water measuring range:	10 % – 76 %
Muscle mass measuring:	13 % – 53 %
Bone measuring range:	2 kg – 10 kg
Size of the product (mm)	310 x 310 x 25

**The manufacturer reserves the right to change the technical specifications and accessories for the respective models.**

**WARNINGS AND SYMBOLS USED ON THE APPLIANCE, PACKAGING OR IN THE INSTRUCTIONS MANUAL:**

HOUSEHOLD USE ONLY. DO NOT IMMERSE IN WATER OR OTHER LIQUIDS.	
	TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, ARRANGES OR PLAYPENS. THIS BAG IS NOT A TOY.
	WARNING
	Please read the instructions manual
	It is not allowed to dispose it together with communal waste
	Do not step onto the scale with wet feet or on wet surface of the scale, there is a risk of slipping.
	AUTO STEP ON – just step on the scale and stand still, the scale will automatically turn ON and start taking measurements.

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© DATE 21/06/2024
Service: Omega electric GmbH, Servicezentrale, Industriering 2, D- 04626 Schmölln.